



IHE INTEGRATION STATEMENT - HYDMEDIA G6 6.02.0

TABLE OF CONTENTS

1	Scope.....	3
2	Purpose.....	3
3	IHE Integration Statement	3
4	Record Revision History	4

1 SCOPE

This document applies to DH Healthcare HYDMedia G6 in Version 6.2.

2 PURPOSE

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below.

3 IHE INTEGRATION STATEMENT

IHE Integration Statement			
Vendor	Product Name	Version	Date
DH Healthcare GmbH	HYDMEDIA G6	6.2	Dec 2020
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
ITI – IT Infrastructure Domain			
Cross-Enterprise Document Sharing (XDS)	Document Repository	None	
	Document Source	Document Replacement	
	Integrated Document Source/Repository	Document Replacement	
	Document Consumer	None	
	Document Administrator	None	
Cross-Enterprise Sharing of Scanned Documents (XDS-SD)	Content Creator	PDF/A	
	Content Consumer	View Option	
Cross-enterprise Document Media Interchange (XDM)	Portable Media Creator	USB, CD-R	
Patient Demographics Query (PDQ)	Patient Demographics Consumer	None	
Document Metadata Subscription (DSUB)	Document Metadata Subscriber	None	
	Document Metadata Notification Recipient	None	
Cross Enterprise User Assertion (XUA)	X-Service User	None	
	X-Service Provider	None	
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	None	

Patient Identifier Cross-referencing HL7 V3 (PIXV3)	Patient Identifier Cross-reference Consumer	None
Consistent Time (CT)	Time Client (TIME_CLIENT)	None
Audit Trail and Node Authentication (ATNA)	Secure Application	None
Sharing Laboratory Reports (XD-LAB)	Content Consumer	None
Cross-Enterprise Sharing of Medical Summaries (XDS-MS)	Content Consumer	None
Exchange of Personal Health Record Content (XPHR)	Content Consumer	View Option
Mobile access to Health Documents (MHD)	Document Source	None
	Document Consumer	None
Remove Metadata and Documents (RMD)	Document Administrator (for HYDMedia XDS Document Source)	None
	Document Repository	None

Internet address for vendor's IHE information: <https://www.dedalusgroup.de/ihe/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.dedalusgroup.de/hl7/
DICOM	https://www.dedalusgroup.de/dicom/

Links to general information on IHE

In North America: www.ihe.net In Europe: www.ihe-europe.net In Japan: www.ihe-j.org/

4 RECORD REVISION HISTORY

Version #	Date	Change description (incl. section)
3	2020-12	Corrected version and document date
1, 2	2020-12	Initial version